

RACK 'EM UP


For safe handling in the
BD RODAC™ Rack



BD RODAC™ racks are specially designed to hold plates tightly in place for convenient carrying in the laboratory and clean room. In addition, you'll have enhanced safety with RODAC™ plate lids held securely in place. RODAC™ racks are available in a variety of colors—blue, orange and green—to differentiate samples.

- Custom made for RODAC™ plates
- Made of welded steel rods and epoxy-coated for maximum protection
- Steam autoclavable at 121°C
- Can be VHP disinfected for isolators

Catalog #	Description
212976	Blue RODAC™ Rack
212977	Orange RODAC™ Rack
212978	Green RODAC™ Rack



BD Diagnostic Systems
7 Loveton Circle
Sparks, MD 21152-0999
800.638.8663
www.bd.com/industrial

BD, BD Logo and RODAC are trademarks of Becton, Dickinson and Company. ©2002 BD.
0-2592 September 2002 Printed in USA